


August 2023
 AC Flow-Vu Medium Mask
 (Yellow Carton)
 P/N 15369-00 ART-0369-00

Mid-York Press Inc
 job# 75341 9/13/23
 Monaghan Medical
 supplied file from artwork
 Small/Medium Mask Cartons
 3+5/8 x 3+3/8 x 7+3/8
 reverse tuck - slit lock
 .016 sbs pgrm# s10589
 approved date

Helps ensure inhaler medication is delivered to the lungs.¹



AeroChamber
Plus **Flow-Vu**
 Anti-Static Valved Holding Chamber
MEDIUM MASK

AeroChamber *Plus* Flow-Vu

Anti-Static Valved Holding Chamber

MEDIUM MASK


For use with Metered Dose or Soft Mist inhalers.



P/N 15369-00 08/2023 ART-0369-00

monaghan

Monaghan Medical Corporation
 Plattsburgh, New York 12901
monaghanmed.com
 518.561.7330 | phone
 800.833.9653 | toll free
 518.561.5088 | fax



COMFORTSEAL
 Mask contours gently to the face providing a secure, comfortable fit.

This device has not been made (or manufactured) with BPA (Bisphenol A) or Natural Rubber Latex.



FLOW-VU Inspiratory Indicator only moves if you have a good mask seal.

For use with Metered Dose or Soft Mist inhalers.


Anti-static Chamber helps maximize aerosol suspension time, maximizing aerosol medication delivery.

Rx ONLY

Easy to use.



Shake the inhaler.



Insert the inhaler into the back of the chamber.



Press inhaler and count for 5 to 6 breaths.



LEARN MORE

Top of the mask should sit on the bridge of the nose and the bottom should rest between the bottom of the lips and chin.

Reference: 1. McIvor RA et al. Can Respir J. 2018 Apr 4;2018:5076259.
 Words or phrases accompanied by "™" and "®" are trademarks and registered trademarks of Monaghan Medical Corporation (MMC) or an affiliate of MMC.
 © 2023 MMC.