


August 2023  
 AC Flow-Vu Small Mask  
 (Orange Carton)  
 P/N 15368-00, ART-0368-00

Mid-York Press Inc  
 job# 75341 9/13/23  
 Monaghan Medical  
 supplied file from artwork  
 Small/Medium Mask Cartons  
 3+5/8 x 3+3/8 x 7+3/8  
 reverse tuck - slit lock  
 .016 sbs pgrm# s10589  
 approved date

Helps ensure inhaler medication is delivered to the lungs.<sup>1</sup>



**AeroChamber**  
*Plus* **Flow-Vu**  
 Anti-Static Valved Holding Chamber  
**SMALL MASK**

**AeroChamber**  
*Plus* **Flow-Vu**

Anti-Static Valved Holding Chamber

**SMALL MASK**


For use with Metered Dose or Soft Mist inhalers.



P/N 15368-00 08/2023 ART-0368-00

**monaghan**

Monaghan Medical Corporation  
 Plattsburgh, New York 12901  
**monaghanmed.com**  
 518.561.7330 | phone  
 800.833.9653 | toll free  
 518.561.5088 | fax



**COMFORTSEAL**  
 Mask contours gently to the face providing a secure, comfortable fit.

This device has not been made (or manufactured) with BPA (Bisphenol A) or Natural Rubber Latex.



**FLOW-VU** Inspiratory Indicator only moves if you have a good mask seal.

For use with Metered Dose or Soft Mist inhalers.


Anti-static Chamber helps maximize aerosol suspension time, maximizing aerosol medication delivery.

**Rx ONLY**

Easy to use.




Shake the inhaler.



Insert the inhaler into the back of the chamber.



Press inhaler and count for 5 to 6 breaths.



LEARN MORE

Top of the mask should sit on the bridge of the nose and the bottom should rest between the bottom of the lips and chin.



Reference: 1. McIvor RA et al. Can Respir J. 2018 Apr 4;2018:5076259.  
 Apr 4, 2018:5076259.  
 Words or phrases accompanied by "™" and "®" are trademarks and registered trademarks of Monaghan Medical Corporation (MMC) or an affiliate of MMC. © 2023 MMC.