



The Monaghan Advantage

- ❖ Breath actuated technology nebulizer
- ❖ Durable, 6-month reusable design
- ❖ Allows for breaks in treatment without drug loss
- ❖ Minimal exposure of caregivers and family members to secondhand aerosols
- ❖ Optimal performance with the OMBRA® Compressor System
- ❖ Nebulizer designed for high level disinfection and dishwasher safe
- ❖ Not made or manufactured with BPA, natural rubber latex, DEHP or lead

AEROECLIPSE® XL R BAN AND OMBRA® TABLE TOP COMPRESSOR

The AEROECLIPSE® XL Reusable Breath Actuated Nebulizer (R BAN) provides a hospital-quality nebulizer treatment for the home care environment. Recommended by clinicians, the AEROECLIPSE® XL R BAN is designed to improve patient compliance for nebulizer treatments outside of the hospital.

Activating only in response to the patient's inspiratory flow, AEROECLIPSE® XL R BAN maintains available drug during exhalation phase, providing patients the ability to take a break or have a conversation without drug loss. Paired with the OMBRA® Compressor, the AEROECLIPSE® XL R BAN puts patients in control of their therapy. Unique and patient friendly features make the system an easy and safe to use, cost effective drug delivery system.

THE FACTS

Recommended by clinicians, the AEROECLIPSE® XL R BAN is designed for nebulizer treatments to improve patient compliance outside of the hospital.

Suitable for use for by all patient populations, the AEROECLIPSE® XL R BAN activates at 15 lpm and can be used with a mouthpiece or our COMFORTSEAL® Masks.

Minimizes exposure of caregivers and family members to dangerous secondhand aerosols, which may carry harmful viruses.


Designed to put patients in control of their therapy at home, AEROECLIPSE® XL R BAN maintains available drug during exhalation phase to give patients time to take a break or have a conversation without drug loss.

MAINTAINS
AVAILABLE
DRUG DURING
EXHALATION
PHASE



518.561.7330 phone
800.833.9653 toll free
518.561.5088 fax

Monaghan Medical Corporation
Plattsburgh, New York 12901
monaghanmed.com



THE FACTS (continued)

Pair the **AEROECLIPSE® XL R BAN** with the **OMBRA® Compressor** for optimal nebulizer performance in a durable home compressor system.* Average treatment time is 6–8 minutes.¹

Reusable for 6-months of single-patient use, the **AEROECLIPSE® XL R BAN** is top shelf dishwasher safe and approved for high-level disinfection.

Kit includes 2 **AEROECLIPSE® XL R BANs** with integrated handle and nebulizer cup holder in a convenient carrying case.

REFERENCE

1. TMI Aerosol Lab, Data on File. 2.5ml, 2.5 mg salbutamol sulphate (Vt = 600 mL, i:E 1:2, 10 bpm) @ 8 minutes

*Warranty information about our **OMBRA®** compressors is available upon request by contacting Monaghan Customer Service at 1-800-833-9653.

CLEANING IS EASY

- **AEROECLIPSE® XL R BAN**, single patient use
 - Designed for 6 months of use at home, clinic or hospital environment
 - Approved for top shelf dishwasher cycle and high level disinfection with bleach, boiling, microwave steam bag or electronic steam sterilizer

ORDERING INFORMATION

PART #	NAME/DESCRIPTION	UNITS PER CASE
69010	AEROECLIPSE® XL Reusable Breath Actuated Nebulizer (R BAN)	cs/10
63506	OMBRA® 120V Table Top Compressor (includes 2 AEROECLIPSE® XL R BANs, 5 Replacement Air Filters, Owner's Manual)	cs/1
66501	OMBRA® Compressor Air Filters, Replacement (5/pack)	cs/1



Contact your Monaghan Sales Territory Manager for more information or a demonstration.

