



The Monaghan Advantage

- ❖ Safe aerosol delivery is possible in breath actuated or continuous modes of nebulizer operation
- ❖ Designed to replace your existing mouthpiece
- ❖ Can be used for adult and pediatric patients age 5 and older
- ❖ Protects caregivers within the treatment environment from breathing in wasted medication or exhaled emissions
- ❖ Can be used in the hospital or home settings

FILTERED MOUTHPIECE KIT

- The only Filter Kit cleared for use with the AEROECLIPSE® BAN™ Nebulizers.
- Filters exhaled emissions, viruses, bacteria and undelivered medication.
- Reduces exposure of potentially harmful emissions to others in the room.
- Direct inhalation pathway ensures no loss of medication for the patient.¹
- Designed to replace the existing mouthpiece.
- Protects caregivers within the treatment environment from breathing in wasted medication or exhaled emissions.
- Tamper resistant valve.
- Standard 22 mm fittings to avoid unnecessary adapters.
- Single patient, disposable device – Mouthpiece approved for 7 days with single use replacement filters.

FILTRATION PERFORMANCE

Viral Filtration Efficiency² >99.99%

Bacterial Filtration Efficiency² >99.999%

REFERENCES

1. Data on file TMI 1872.

2. The bacterial and viral filtration efficiency performance of the filtration media has been tested by the material's manufacturer. Report available upon request.

ORDERING INFORMATION

PART #	NAME/DESCRIPTION	UNITS PER CASE
60050	FILTERED MOUTHPIECE KIT Filter Kits include a one-way valved tee with mouthpiece and filter assembly.	cs/50
60125	AEROECLIPSE® II BAN™ with Universal and Filtered Mouthpiece Kit	cs/25
60225	AEROECLIPSE® II BAN™ with EZ Twist Tubing and Filtered Mouthpiece Kit	cs/25
60100	Filter Alone, Replacement	cs/100

Contact your Monaghan Sales Territory Manager for more information or a demonstration.



REGISTERED
ISO 13485
QUALITY SYSTEM